

Reprocessed for



Sustainability Solutions

Instructions for Use Reprocessed Posey Bed Alarm Sensor Pad

Reprocessed Device for Single Use. Reprocessed for Stryker Sustainability Solutions.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

NON-STERILE NOT MADE WITH NATURAL RUBBER LATEX

Explanation of Icons



Date of Reprocessing



Use by Date

REF

Product Code



Do Not Reuse



See Instructions For Use

Indications for Use

The Stryker Sustainability Solutions Reprocessed Posey Bed Alarm Sensor Pads alert staff that a patient has taken weight off the sensor and may be at immediate risk of falling. However, fall alarm systems do not prevent falls and are not a substitute for good nursing care and regular visual monitoring.

These sensor pads must be connected to a compatible Posey Bed Alarm Unit.

Instructions for Use:

1. Position sensor pad on top of mattress, across width of bed, in one or two positions (see Figure):
 - a. Directly under patient's shoulders; alarm will sound when patient sits up.
 - b. Directly under patient's buttocks; alarm will sound when patient moves off the sensor as in bed exit.
2. Secure the sensor pad to the mattress using the metal clips.
3. Place the linens on the bed over the sensor pad.
4. Connect the sensor pad to the Posey alarm unit.

System Test:

1. Plug sensor into alarm unit.
2. Turn alarm unit on and verify power light is lit.
3. Check to make sure the audio alarm is functioning properly by lifting weight off the sensor pad and activating the alarm.
4. Ensure the nurse call interface, when used, is connected properly and is activated properly.

Reprocessed Posey Bed Alarm Sensor Pad

Warnings

To avoid possible injury or death:

- Read sensor pad and alarm instructions before use.
- This device may not be suitable for use with all “high fall risk” patients. Restraints may still be required.
- This device should not be a substitute for maintaining the routine visual monitoring protocol by healthcare personnel.
- Never connect the sensor pad to other manufacturers’ alarm units.
- Always test system per instructions prior to use.
- Be sure restraint straps do not cross over the sensor pad! Restraint straps may prevent alarm activation if patient falls off mattress.
- When using a restraint, always use the side rails in the “up” position. The patient’s body or limb should not fit over, under, around, through or between side rails.
- The alarm may fail to sound if the pad is bent or folded.

Always verify system is working properly before leaving the patient unattended.

Do not submerge the sensor pad in liquid or sterilize it with heat.

Cleaning: wipe clean with liquid disinfectant.

Do not remove labeling from sensor pad during the life of the product.

Do not use with memory foam or low air-loss mattresses or overlays. A foam pad on top of the mattress may diffuse patient weight and prevent sensor from activating. If alarm does not function, place sensor pad above foam pad and re-test.




Figure – 1 Sensor Pad Locations

The OEM information listed on the label is provided as device ID prior to reprocessing and may contain the trademarks of unrelated third parties that do not sponsor this device.

This product and its packaging have been exposed to ethylene oxide gas (EtO). Even though the product is processed in compliance with all applicable laws and regulations relating to EtO exposure, proposition 65, a State of California voter initiative, requires the following notice:

Warning: this product and its packaging have been exposed to ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.