# IntelliVue Reusable ECG Lead Sets and Trunk Cables

## Instructions for Use

3 ECG	4 ECG	5 ECG	5 (3+2) ECG	6 (4+2) ECG	ECG
AAMI IEC M1671A, 75A* M1672A, 78A*	AAMI IEC M1532A, 57A* M1533A, 58A*	AAMI IEC M1973A*, 76A, 79A* M1978A, 84A*	AAMI IEC M1968A M1971A, 74A*	AAMI IEC M1680A, 84A* M1681A, 85A*	
AAMI IEC M1673A M1674A	AAMI IEC M1537A M1538A	AAMI IEC M1602A M1604A	AAMI IEC M1644A M1645A	AAMI IEC M1682A M1683A	
AAMI IEC M1622A, 24A M1626A	AAMI IEC 989803192241 989803192251	AAMI IEC 989803192181 989803192191	AAMI IEC M1647A M1648A	AAMI IEC 989803192221 989803192231	
			AAMI IEC 989803192201 989803192211		
					M1663A (5+5) M1665A (6+4) M1667A (6) M1668A (5) M1669A (3) M1949A (5+5) 989803170171 (3)* 989803170181 (5)*

\* Indicates products approved for use in Operating Room (OR) environments.

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

**[www.philips.com/healthcare/eIFU](http://www.philips.com/healthcare/eIFU)**


(BG) Указания за употреба .....	
(CS) Návod k použití .....	
(DA) Brugervejledning .....	
(DE) Gebrauchsanweisung .....	
(EL) Οδηγίες χρήσης .....	
(ET) Kasutusjuhend .....	
(FI) Käyttöohjeet .....	
(HU) Használati útmutató .....	
(IT) Istruzioni d'uso .....	
(LT) Naudojimo instrukcija .....	
(LV) Lietošanas instrukcija .....	
(NL) Gebruiksaanwijzing .....	
(NO) Bruksanvisning .....	
(PL) Instrukcja obsługi .....	
(RO) Instrucţiuni de utilizare .....	
(SK) Návod na obsluhu .....	
(SL) Navodila za uporabo .....	
(SV) Bruksanvisning .....	
(TR) Kullanım Talimatları .....	

AAMI	3 ECG	4 ECG	5 ECG	6 ECG				OR		REF
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>			M1671A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1675A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>			M1673A
<input checked="" type="checkbox"/>					0.7 m		<input checked="" type="checkbox"/>			M1624A
<input checked="" type="checkbox"/>					0.45 m		<input checked="" type="checkbox"/>			M1622A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1532A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		M1557A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1537A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192241
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>			M1968A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1976A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		M1979A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1973A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1602A
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>			M1644A
<input checked="" type="checkbox"/>					0.7/1.3 m		<input checked="" type="checkbox"/>			M1647A
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192201
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	989803192181
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			M1680A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1684A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			M1682A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192221

IEC	3 ECG	4 ECG	5 ECG	6 ECG				OR		REF
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>			M1672A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1678A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>			M1674A
<input checked="" type="checkbox"/>					0.7 m		<input checked="" type="checkbox"/>			M1626A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1533A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		M1558A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1538A
<input checked="" type="checkbox"/>					1.0 m		<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192251
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>			M1971A
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1974A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1978A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		M1984A
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>			M1645A
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>				M1604A
<input checked="" type="checkbox"/>					0.7/1.3 m		<input checked="" type="checkbox"/>			M1648A
<input checked="" type="checkbox"/>					1.0/1.6 m		<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192211
<input checked="" type="checkbox"/>					1.0 m	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>	989803192191
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			M1681A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		M1685A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			M1683A
<input checked="" type="checkbox"/>					1.0/1.6 m	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	989803192231

AAMI + IEC		OR	REF
AAMI + IEC (3)	2.7 m		M1669A
AAMI + IEC (5)	2.7 m		M1668A
AAMI + IEC (6)	2.7 m		M1667A
AAMI + IEC (5+5)	2.0 m		M1663A
AAMI + IEC (5+5)	2.7 m		M1949A
AAMI + IEC (6+4)	2.7 m		M1665A
AAMI + IEC (3)	2.7 m	<input checked="" type="checkbox"/>	989803170171
AAMI + IEC (5)	2.7 m	<input checked="" type="checkbox"/>	989803170181

REF 989803193821

## INDICATIONS FOR USE

Philips IntelliVue reusable ECG lead sets and trunk cables are limited by the indications for use of the connected monitoring and diagnostic equipment, and are for use only by healthcare professionals. They are designed for multiple patient use and indicated for monitoring cardiac signals for both diagnostic and monitoring purposes in adult, pediatric, and neonatal patients.

Philips operating room (OR) ECG trunk cables are indicated for use to monitor patient ECG in an electrosurgery (ESU) environment when used in combination with a compatible ECG patient lead set on adult, pediatric, infant, and neonatal ECG applications.

## DESCRIPTION

Tables on pages 1- 2 define physical characteristics and applications for all ECG lead set and trunk cable products included in this IFU. A descriptive key to all column headings within those tables is provided below.

<b>AAMI</b>	Check mark ✓ indicates that the lead set is AAMI color coded.
<b>IEC</b>	Check mark ✓ indicates that the lead set is IEC color coded.
<b>3 ECG</b> <b>4 ECG</b> <b>5 ECG</b> <b>6 ECG</b>	Indicates the total number of ECG leads in each lead set, and the type of electrode connector at the end of each lead:  = Grabber  = Snap  = Mini Clip  = DIN (each DIN lead requires an electrode connector adapter [REF 989803193821]).
	Identifies length(s) of lead set leads, or trunk cable.
	Check mark ✓ identifies lead sets that are appropriate for chest application.
	Check mark ✓ identifies lead sets that are appropriate for arm/leg application.
	Check mark ✓ identifies lead sets and/or trunk cables that may be used in operating room environments.
	Check mark ✓ identifies lead sets that require electrode connector adapters. One adapter ([REF 989803193821]) is needed for each DIN lead. See WARNINGS section.
	Defines the lead set or trunk cable catalog reference number.
<b>AAMI + IEC (x) or (x-x)</b>	Defines ECG trunk cables as being both AAMI and IEC compatible, with (x) or (x-x) identifying each trunk cable's lead set connector pin configuration.

## COMPATIBILITY

IntelliVue ECG lead sets and trunk cables can be used with any ECG monitor or defibrillator for which they are listed as accessories in that product's IFU.

## DISPOSAL

Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

## WARNINGS

- Leads that have DIN  electrode connectors require a separate electrode connector adapter for each lead. For more information reference the *Reusable Electrode Adapter Instructions for Use*.
- Refer to the IFU for the monitor/defibrillator for additional warnings and cautions.
- Refer to the IFU for the monitor/defibrillator for information regarding proper lead/electrode placement that complies with standard AAMI or IEC practices.
- Ensure that the patient is properly grounded during ESU procedures.
- Do not use ECG leads and cables in magnetic resonance imaging (MRI) environments or during MRI procedures.
- Ensure that ECG leads and cables are carefully positioned to avoid entanglement, choking, and strangulation.
- When used to monitor ECG in the OR, ensure that you use accessories that provide appropriate ESU protection (orange OR lead sets, or ICU lead sets paired with Philips OR trunk cables 989803170171 or 989803170181).
- OR lead sets and/or OR trunk cables cannot be used to measure respiration.

## CAUTIONS

- Do not use if visual inspection reveals damage or contamination.
- Do not use lead sets/cables in excessively wet environments.
- Ensure that the trunk cable is properly plugged into the monitor/defibrillator and lead set.

## CLEANING AND DISINFECTION

### Lead Set and Trunk Cable Cleaning

To reduce risk of cross-infection, always clean and disinfect the product before applying it to a different patient.

Modestly dampen a lint-free cloth with one of the approved cleaning/disinfecting agents listed in the table that follows and gently wipe all cable and connector surfaces. After cleaning, wipe the outside surfaces of the cable/leads/connectors with a clean lint-free cloth dampened with clean water.

**Note:** Do not immerse the trunk cable connectors or lead set connectors in any fluid.

**Note:** Follow the manufacturer's instructions regarding duration of cleaner application.

### DIN Electrode Connector Adapter Cleaning

Approved cleaning and disinfecting agents and instructions for cleaning the DIN electrode connector adapter are provided in the *Reusable Electrode Adapter Instructions for Use*.

## STERILIZATION

### Lead Set and Trunk Cable Sterilization

Use Ethylene Oxide (EtO) sterilization only. Follow all safety precautions regarding aeration after EtO exposure.

### DIN Electrode Connector Adapter Sterilization

Reference the *Reusable Electrode Adapter Instructions for Use*.

Lead Set and Trunk Cable Cleaning and Disinfecting Agents (see Notes below table)	
Isopropanol	☑
Bacillo <sup>®</sup> AF	☑
Bacillo <sup>®</sup> 25	☑
Melisepto <sup>®</sup>	☑
Accel TB RTU	☑
Oxivir <sup>®</sup> Tb Cleaner Disinfectant	☑
Oxivir <sup>®</sup> Tb Wipes	☑
Carpe Diem <sup>TM/MC</sup> Tb	☑
Super Sani-Cloth	☑
SANI-CLOTH <sup>®</sup> PLUS Germicidal Disposable Wipes	☑
SANI-CLOTH <sup>®</sup> HB Germicidal Disposable Wipes	☑
Hexanios G+R	☒
Sodium hypochlorite, bleach	☒
Surfanios Citro	☒

☑ = Approved cleaner/disinfectant.

☒ = Cleaning/disinfecting agent is known to cause damage to plastic material and should be avoided, if possible.

After cleaning/disinfecting, ensure that all residue is removed by wiping the surface using a clean lint-free cloth dampened with clean water.

Trademark indicator <sup>®</sup> or <sup>TM/MC</sup> defines a trademarked product name of the respective manufacturer.

## SYMBOL DEFINITIONS

Symbols that may appear on products or product packaging are defined below.

= Quantity per pouch

Refer to Instruction Manual/Booklet

Attention, consult accompanying documents

Consult Instructions for Use

Defibrillation proof type CF Applied Part

**R<sub>x</sub> only**

Prescription use only

Not manufactured with natural rubber latex

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To locate your local Philips sales office, visit [www.healthcare.philips.com](http://www.healthcare.philips.com) or contact a Philips regional office listed below.

**United States:**

Philips Healthcare  
North America Corporation  
3000 Minuteman Road  
Andover, MA 01810-1099  
+1 800 225-0230 toll free

**Canada:**

Philips Healthcare Canada  
281 Hillmount Road  
Markham, ON L6C 2S3  
(800) 291-6743

**Latin America Headquarters:**

Philips Healthcare Brazil  
Rua Verbo Divino, 1400 - 7º andar  
Chacara Santo Antonio  
São Paulo, SP 04719-002 Brazil  
Tel: +55 11 2125 0744 Fax: +55 11 2125 0761

**Europe, Middle East, and Africa:**

Philips Medizin Systeme  
Böblingen GmbH  
Hewlett-Packard Str. 2  
71034 Böblingen, Germany  
Tel: (+49) 7031 463 2254 Fax: (+49) 7031 463 1552

**Asia Pacific Headquarters:**

Philips Electronics Singapore Private Limited  
Philips Healthcare  
620A Lorong 1 Toa Payoh  
Singapore 319762  
Tel: 1-800-PHILIPS Fax: +65 6255 4853

 989803170171, 70181, 92181,  
92191, 92211, 92201, 92221, 92231,  
92241, 92251

Philips Medical Systems  
3000 Minuteman Road  
Andover, MA 01810-1099, USA

All other lead set and trunk cable  
product part numbers in this document  
that are not listed above.

Philips Medizin Systeme  
Böblingen GmbH  
Hewlett Packard Str. 2  
71034 Böblingen, Germany

For more information or to reorder, go to  
[www.philips.com/healthcarestore](http://www.philips.com/healthcarestore).

Subject to modification

Copyright © 2014  
Koninklijke Philips N.V.  
September 2014 Edition

All rights reserved  
IFU Part No. 453564484201

*A-453564484201-2 Rev. A*