

Reprocessed by



HEALTH SERVICES

434 Industrial Lane Birmingham, AL 35211  
1-866-943-6670

## **Instructions for Use**

### **Hygia Reprocessed Aircast Elite VenaFlow Calf & Thigh Sleeve IFU**

#### **Reprocessed Device for Single Use.**

**Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

**NON-STERILE      NOT MADE WITH NATURAL RUBBER LATEX**

#### **Sleeves Description**

Hygia Reprocessed Aircast Elite Venaflow Calf & Thigh Sleeves

#### **Indications for Use**

The Hygia Health Services Reprocessed Aircast Venaflow Elite sleeves are designed to increase venous blood flow in at-risk patients in order to prevent deep vein thrombosis. The sleeves are a non-invasive therapeutic method for prevention of deep vein thrombosis. The sleeves are also used in the treatment of venous leg/foot ulcers and edema which are disorders associated with venous insufficiency.

#### **Warnings/Precautions**

- Proper sleeve positioning must be assured to avoid the potential for pressure points on the leg.
- Ensure proper connections to the Venaflow pump and sleeves.
- Kinked or twisted tubing may restrict air flow.
- Sleeves should be removed if patient experiences numbness, tingling, or leg pain.
- To minimize local air movement when using sleeves in operating room, cuff cooling should be turned off, if available.

#### **Directions for Use**

1. Remove Aircast sleeves from plastic bag.
2. Apply each sleeve with the aircell centered on back of the leg and with the tubes toward the foot.
3. Secure sleeve straps snug, but not tight.
4. Be sure to check that all tube assembly connections are tight and that sleeves are snug to achieve full pressure. Always turn off pump before readjusting or removing device.
5. Refer to original equipment manufacturer manual for further instructions.

#### **Contraindications**

1. Any Local leg condition in which sleeves/cuffs would interfere such as dermatitis, vein ligation (immediate postoperative), gangrene, or recent skin graft
2. Severe arteriosclerosis or other ischemic vascular disease
3. Massive edema of legs or pulmonary edema from congestive heart failure
4. Extreme deformity of leg
5. Suspected existing deep venous thrombosis
6. Congestive heart failure

The OEM information listed on the label is provided as device ID prior to reprocessing and may contain the trademarks of unrelated third parties that do not sponsor this device.

This product and its packaging have been exposed to ethylene oxide gas (EtO). Even though the product is processed in compliance with all applicable laws and regulations relating to EtO exposure, proposition 65, a State of California voter initiative, requires the following notice:

Warning: this product and its packaging have been exposed to ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.