

Reprocessed by



434 Industrial Lane Birmingham, AL 35211
1-866-943-6670

Instructions for Use Hygia Stethoscope IFU

Reprocessed Device for Single Use.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

NON-STERILE NOT MADE WITH NATURAL RUBBER LATEX

Indications:

The Hygia Health Services Reprocessed Single Head and Dual Head Stethoscope is a non-invasive device used to transmit to the examiner's ear sounds produced in the body.

Instructions for Use:

1. Inspect the stethoscope for quality and damage.
2. Find a relatively quiet area to ensure that the body sounds are not overpowered by background noise.
3. Insert each of the two ear pieces into your ears. Be sure that the ear pieces fit snugly and have a good seal to keep out ambient noise.
4. Hold the chestpiece (diaphragm or bell) directly against the skin for best sound transmission.

Warnings:

- Avoid extreme heat, cold, solvents and oils.
- Do not leave exposed to direct sunlight for prolonged periods of time.
- The entire stethoscope including the chestpiece, diaphragm, tubing and eartips can be cleaned with isopropyl alcohol or soapy water.
- After cleaning, make sure to thoroughly rinse and dry all parts prior to storage or use.
- Never immerse the stethoscope in any liquid or subject it to any sterilization methods

The OEM information listed on the label is provided as device ID prior to reprocessing and may contain the trademarks of unrelated third parties that do not sponsor this device.

This product and its packaging have been exposed to ethylene oxide gas (EtO). Even though the product is processed in compliance with all applicable laws and regulations relating to EtO exposure, proposition 65, a State of California voter initiative, requires the following notice:

Warning: this product and its packaging have been exposed to ethylene oxide. The packaging may expose you to ethylene oxide, a chemical known to the State of California to cause cancer or birth defects or other reproductive harm.